

Index

Sr. No.	Experiments	Page No.	Date	Remark	Signature
Quantitative estimation (Long experiment)					
1	To determine the concentration in terms of molarity of KMnO_4 by titrating it against (0.1 M) standard solution of oxalic acid.	5 TO 7	5/10/23		
2	To determine the concentration in terms of molarity of KMnO_4 by titrating it against (0.1 M) standard solution of ferrous ammonium sulphate.	9-11	5/10/23		
Chemical kinetics (Long experiment) (Any two)					
3	To study the effect of concentration on rate of reaction between $\text{Na}_2\text{S}_2\text{O}_3$ and HCl.	13-17	26/10/23		
4	To study the effect of temperature on rate of reaction between $\text{Na}_2\text{S}_2\text{O}_3$ and HCl.				
5	To study the rate of reaction between KIO_3 and $\text{Na}_2\text{S}_2\text{O}_3$ using starch as an indicator.				
Electrochemistry					
6	Study of variation of cell potential of cell $\text{Zn} \text{Zn}^{2+} \text{Cu}^{2+} \text{Cu}$ with change in concentration of electrolyte CuSO_4 and ZnSO_4 at room temperature.				
Thermochemistry (Short experiment) (Any three)					
7	To determine enthalpy of dissolution of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ in water at room temperature.				
8	To determine enthalpy of neutralization of HCl and NaOH.	35 TO 37	6/10/23		
9	To determine enthalpy change of displacement of Cu from CuSO_4 by Zn.	39 TO 41	6/10/23		
10	To determine enthalpy change during the interaction between acetone and chloroform				
Inorganic preparations (Short experiment) (Any two)					
11	To prepare a pure sample of ferrous ammonium sulphate.				
12	To prepare a pure sample of potash alum.				
13	To prepare a pure sample of potassium trioxalatoferate (III) complex.				

Sr. No.	Experiments	Page No.	Date	Remark	Signature
Organic preparations (Short experiment) (Any two)					
14	To prepare dibenzalacetone.				
15	To prepare 2-naphthol aniline dye.				
16	To prepare acetanilide from acetyl chloride.				
17	To prepare p-nitro acetanilide from acetanilide.				
Identification of functional group (Any four compounds)					
18	Test for functional group present in given organic compound.	71	27/10/23		
Characteristic tests of carbohydrates, proteins, fats and oil in pure samples and their detection in given food stuffs. (Any one of each)					
19	Characteristic tests of carbohydrate, proteins, fats and oils.	77	27/10/23		
Qualitative analysis for two acidic radicals (Any three water soluble salt mixture)					
Detection of two acidic radicals qualitatively from given inorganic salt mixture.					
20	The given inorganic salt mixture no.1 contains following two acidic radicals. i CO_3^{2-} (Carbonate) ii Cl^- (Chloride)	93 TO 94.	12/10/23		
21	The given inorganic salt mixture no.2 contains following two acidic radicals. i Cl^- (Chloride) ii SO_4^{2-} (Sulphate)	95 TO 96.	12/10/23		
22	The given inorganic salt mixture no.3 contains following two acidic radicals. i CO_3^{2-} (Carbonate) ii Br^- (Bromide)	97 TO 98.	27/10/23		
Qualitative analysis for two basic radicals (Any three water soluble salt mixture)					
Detection of two basic radicals qualitatively from given inorganic salt mixture					
23	The given inorganic salt mixture no.1 contains following two basic radicals. i NH_4^+ (Ammonium) ii Cu^{2+} (Copper ion)	113-113	13/10/23		
24	The given inorganic salt mixture no.2 contains following two basic radicals. i NH_4^+ (Ammonium) ii Ni^{2+} (Nickel)	114-116	19/10/23		
25	The given inorganic salt mixture no.3 contains following two basic radicals. i NH_4^+ (Ammonium) ii Zn^{2+} (Zinc)	117-119	19/10/23		

Sr. No.	Experiments	Page No.	Date	Remark	Signature
List of Activities (Any four)					
1	Preparation of phthalic anhydride from phthalic acid.				
2	Preparation of succinic anhydride from succinic acid.				
3	Preparation of methyl salicylate.				
4	Measurement of e.m.f. of commercial cells by using voltmeter. (Any four cells)				
5	Identification of Resin identification code (RIC) on given different plastic materials. (Any two plastic material)	128-129	5/10/23		
6	Preparation of 100 mL 0.1 M standard solution of ferrous ammonium sulphate. (Mohr's salt)	130-131	5/10/23		

Important Instructions

- 1) Read instructions of do's and don'ts carefully.
- 2) Record your observations and calculations on space provided in the notebook.
- 3) Take signature of the batch incharge after completion of experiment on same day.
- 4) Read carefully the begining of each section, a concise explanatory account of theoretical experiments has been given.
- 5) Solve exercise of each experiment after completion of experiment on same day.
- 6) Use given space provided after each experiment for notes, rough work and rough calculations and then fair it to main experiment.
- 7) A Scheme of marking for practical examination is given on **page no. 143** for accurate working.
- 8) A format of question paper for annual practical examination is given for your information on **page no. 145**.
- 9) At the end, tables of logarithms is given for calculation.
- 10) At the end, tables for preparation of different concentration solutions and reagents are given in **appendix I to V**.


Observations:

1. Solution in a burette : KMnO_4
2. Solution by a pipette : oxalic acid
3. Solution in conical flask : 10mL oxalic acid + dil. H_2SO_4 (1 test tube)
4. Indicator : KMnO_4 acid
5. End point : Colourless to light pink
6. Chemical Equation : $2\text{KMnO}_4 + 5\text{H}_2\text{C}_2\text{O}_4 + 3\text{H}_2\text{SO}_4 \rightarrow \text{K}_2\text{SO}_4 + 2\text{MnSO}_4 + 8\text{H}_2\text{O} + 10\text{CO}_2 \uparrow$

Observation Table:

Burette level	Pilot reading	Burette reading in mL			C.B.R.
		I	II	III	
Final	9	9.5	9.5	9.5	9.5 (x) mL
Initial	to 10	0.0	0.0	0.0	
Difference	mL	9.5	9.5	9.5	

Calculation: From above chemical equation



$$2 \times 158 \text{ g (316 g)} \equiv 1000 \text{ mL } 5 \text{ M} \text{ (molar mass of } \text{KMnO}_4 = 158 \text{ g/mol)}$$

$$\therefore 1000 \text{ mL } 5 \text{ M oxalic acid} \equiv 316 \text{ g of } \text{KMnO}_4$$

$$\therefore 10 \text{ mL } 0.1 \text{ M oxalic acid} = \frac{316 \times 10 \times 0.1}{1000 \times 5}$$

$$= 0.0632 \text{ g of } \text{KMnO}_4$$

\therefore Hence 9.5 (x) C.B.R. mL of KMnO_4 solution contains = 0.0632 g of KMnO_4

$$\therefore 1000 \text{ mL of } \text{KMnO}_4 \text{ contains} = \frac{0.0632 \times 1000}{9.5(x) \text{ CBR}}$$

$$\text{Hence molarity of } \text{KMnO}_4 \text{ in the solution is} = \frac{0.0632 \times 1000}{9.5(x) \text{ CBR} \times 158}$$

$$= \frac{0.4}{9.5(x) \text{ CBR}}$$

$$= 0.0421 \text{ M}$$

Space for log calculation

Number	log
9.5	0.977

Result: The molarity of KMnO_4 solution is 0.0421 M

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. The molar mass of oxalic acid is
 - a. 126 u
 - ✓ b. 126 g/mol
 - c. 12.6 u
 - d. 12.6 g/mol
2. To prepare 1 M solution of oxalic acid
 - ✓ a. 126 g of oxalic acid dissolve in distilled water and diluted to 1 L
 - b. 12.6 g of oxalic acid dissolve in distilled water and diluted to 1 L
 - c. 1.26 g of oxalic acid dissolve in distilled water and diluted to 1 L
 - d. 0.126 g of oxalic acid dissolve in distilled water and diluted to 1 L
3. Water of crystallization present in oxalic acid is/are.....
 - a. 1
 - ✓ b. 2
 - c. 3
 - d. 5
4. Oxalic acid is used to prepare standard solution because it is a.....
 - a. substance
 - ✓ b. primary standard substance
 - c. secondary standard substance
 - d. tertiary standard substance
5. The quantity of oxalic acid required to prepare 0.1 M 100 mL standard solution of oxalic acid is
 - a. 126 g
 - b. 12.6 g
 - ✓ c. 1.26
 - d. 0.126 g

Short answer questions

1. Calculate the molar mass of oxalic acid. ($H = 1, C = 12, O = 16$)

Ans.....

Molar mass of oxalic acid is 126 g/mol.

2. Why heating is required in oxalic acid and potassium permanganate titration?

Ans.....

The heating is required in oxalic acid & potassium permanganate titration because reaction is very slow at room temperature.

3. Why one test tube of 2M sulphuric acid is required in permanganate titration?

Ans.....

One test of dilute H_2SO_4 is required in permanganate solution because it prevent oxidation of manganese dioxide.

4. What is the oxidation state of carbon atom in oxalic acid after completion of redox titration?

Ans.....

The oxidation state of carbon atom in oxalic acid after completion redox titration is +4.

5. Why primary standard substances are used to make standard solutions?

Ans.....

Primary standard is a substance which is very pure representative of no. of moles of substance & easily weighted used to determine as unknown concn of $KMnO_4$.

Remark and sign of teacher:


Observations:

1. Solution in a burette : KMnO₄
2. Solution by a pipette : 0.1 M F.A.S.
3. Solution in conical flask : 10 mL 0.1M F.A.S + 1 test tube dil. H₂SO₄
4. Indicator : KMnO₄
5. End point : Colourless to light pink
6. Chemical Equation : $2\text{KMnO}_4 + 10\text{FeSO}_4(\text{NH}_4)_2\text{SO}_4 + 8\text{H}_2\text{SO}_4 \rightarrow \text{K}_2\text{SO}_4 + 2\text{MnSO}_4 + 5\text{Fe}_2(\text{SO}_4)_3 + 10(\text{NH}_4)_2\text{SO}_4 + 8\text{H}_2\text{O}$

Observation Table:

Burette level	Pilot reading	Burette reading in mL			C.B.R.
		I	II	III	
Final	10	10.5	10.5	10.5	10.5 (x) mL
Initial	to	0.0	0.0	0.0	
Difference	mL	10.5	10.5	10.5	

Calculation: From above chemical equation



$2 \times 158 = 10$ moles of mohr's salt (molar mass of KMnO₄ = 158 g)

$\therefore 1000 \text{ mL } 10\text{M Mohr's salt} \equiv 316 \text{ g of KMnO}_4$

$$\therefore 10 \text{ mL } 0.1\text{M Mohr's salt} = \frac{316 \times 10 \times 0.1}{1000 \times 10}$$

$$= 0.0316 \text{ g of KMnO}_4$$

\therefore Hence 10.5 (x) CBR mL of KMnO₄ solution contains = 0.0316 g of KMnO₄

$$\therefore 1000 \text{ mL of KMnO}_4 \text{ contains} = \frac{0.0316 \times 1000}{10.5(x) \text{ CBR} \times 158}$$

$$\begin{aligned} \text{Hence molarity of KMnO}_4 \text{ in the solution is} &= \frac{0.0316 \times 1000}{10.5(x) \text{ CBR} \times 158} \\ &= \frac{0.2}{10.5(x) \text{ CBR}} \\ &= 0.019 \text{ M} \end{aligned}$$

Result: The molarity of KMnO₄ solution is 0.019 M

Space for log calculation

Number	log
10.5	

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. The indicator used in Redox titration is----
a. phenolphthalein b. methyl orange
c. methylene blue ✓d. potassium permanganate
2. Select the titrant used in redox titration
a. oxalic acid b. F.A.S ✓c. KMnO_4 d. H_2SO_4
3. The oxidation state of Mn in KMnO_4 in redox titration before titration is +7, the oxidation state of Mn after completion of titration is----
✓a. +2 b. +3 c. +4 d. +7
4. The oxidation state of Iron in Mohr's salt is +2, the oxidation state of Iron after completion of titration is----
a. +1 b. +2 ✓c. +3 d. +4
5. The role of Mn^{2+} ion during oxidation-reduction titration of KMnO_4 by oxalic acid is----
✓a. catalyst b. reductant c. oxidant d. reactant

Short answer questions

1. What specific name is given to the permanganate titrations?

Ans..... Redox titration is given to the permanganate solutions.

2. Why heating is not required in F.A.S. and potassium permanganate titration?

Ans..... In this titration heating is not required in F.A.S & Potassium permanganate titration because the reaction rate is very high at Room temperature.

3. Why is dilute sulphuric acid added while preparing a standard solution of Mohr's salt?

Ans..... Simplest reason for using dil. H_2SO_4 added while preparing a standard solution of Mohr's salt to get clear solution & prevent hydrolysis of ferrous Sulphate.

4. Calculate the molar mass of Mohr's salt ($\text{Fe} = 56, \text{S} = 32, \text{N} = 14, \text{H} = 1, \text{O} = 16$)

Ans..... Mohr's salt ($\text{FeSO}_4 \cdot (\text{NH}_4)_2\text{SO}_4 \cdot 6\text{H}_2\text{O}$)

$$= 56 + 32 + 64 + 28 + 8 + 32 + 64 + 12 + 96$$

$$= 392 \text{ g/mol}$$

5. Calculate the amount of F.A.S. required to prepare M/20 standard solution of F.A.S.

Ans..... M/20 Standard of F.A.S. means 0.05 Solt.

$$1 \text{ M FAS} = 392 \text{ gm}$$

$$0.1 \text{ M FAS} = 39.2 \text{ gm}$$

$$0.05 \text{ M FAS} = 19.6 \text{ gm.}$$

Remark and sign of teacher:

Procedure:

1. Take four conical flasks (250 mL), wash with water and label them as A, B, C, D respectively.
2. Using 50 mL burette take exactly 20, 30, 40 and 50 mL 0.1 M $\text{Na}_2\text{S}_2\text{O}_3$ in the flasks A, B, C, D respectively.
3. Using another 50 mL burette add 28, 18 and 8 distilled water to flasks A, B and C respectively. (There is no addition of distilled water in flask D)
4. With the help of 10 mL measuring cylinder, add 2 mL 1 M HCl to the flask A containing 20 mL 0.1 M $\text{Na}_2\text{S}_2\text{O}_3$ and 28 mL distilled water. Start the stop watch immediately.
5. Shake and keep the conical flask on a paper having cross mark and view the cross mark through reaction mixture from top of the conical flask.
6. When the cross mark on the paper just become invisible stop the stop watch immediately and record the time required in seconds.
7. Repeat the experiment by adding 2 mL 1 M HCl to flasks B, C and D, simultaneously. By using stop watch record the time required, when cross mark on the paper just becomes invisible. (Use same paper having cross mark for flasks B, C and D)

Observation Table:

Conical Flasks	0.1 M $\text{Na}_2\text{S}_2\text{O}_3$ in mL	H ₂ O mL	1M HCl mL	Time required to become cross mark invisible (t) sec	1/t s ⁻¹	Concentration of $\text{Na}_2\text{S}_2\text{O}_3$
A	20	28	2	414(s)	0.0024	0.04
B	30	18	2	233(s)	0.0043	0.06
C	40	8	2	146(s)	0.0068	0.08
D	50	-	2	138(s)	0.0072	0.1

Graph: Plot a graph of 1/t against concentration of sodium thiosulphate solution.

Result:


Rate of reaction (1/t) is directly proportional to conc. of $\text{Na}_2\text{S}_2\text{O}_3$

Remark and sign of teacher:

Scale :

$$y \text{ axis} = 1 \text{ cm} = 0.001 \text{ s}^{-1}$$

$$x \text{ axis} = 1 \text{ cm} = 0.01 \text{ cm} = \\ 0.01 \text{ mol/dm}^3$$



MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. The rate of reaction between 0.1 M $\text{Na}_2\text{S}_2\text{O}_3$ and 1 M HCl does NOT depends on
 - a. temperature
 - b. concentration
 - c. pressure
 - d. catalyst
2. Sodium thiosulphate solution reacts with hydrochloric acid solution to produce
 - a. a colloidal solution of sulphur
 - b. transparent solution of sulphur
 - c. black solution of sulphur
 - d. white precipitate of NaCl
3. The time required for completion of the reaction, when we add 2 mL of 1M HCl to the flasks A, B, C and D containing increasing order of concentration of 0.1 M $\text{Na}_2\text{S}_2\text{O}_3$
 - a. increases
 - b. decreases
 - c. first decreases and then increases
 - d. first increases and then decreases
4. The nature of graph for effect of concentration on rate of reaction between 0.1 M $\text{Na}_2\text{S}_2\text{O}_3$, and 1 M HCl is a straight line with
 - a. decreasing slope
 - b. increasing slope
 - c. increasing slope intersecting to y-axis
 - d. decreasing slope intersecting to x-axis

Short answer questions

1. Define rate of reaction.

Ans. Rate of reaction can be measured in terms of either decrease in conc. of any one of the reactants or increase in concn of any one of the products per unit time.

2. Mention the factors affecting the rate of reaction.

Ans. Factors affecting the rate of reacn is temperature and catalyst.

3. Write the name of law for the study of effect of concentration on the rate of reaction.

Ans. A differential rate law express reaction rate Inter. of changes in the concn of one or more reactants. Name of law — law of mass action.

4. Explain the nature of graph of $1/t$ against concentration of $\text{Na}_2\text{S}_2\text{O}_3$ solution.

Ans. Rate of Reacn is directly proportional to concn of $\text{Na}_2\text{S}_2\text{O}_3$. The nature of graph $1/t$ against concn of $\text{Na}_2\text{S}_2\text{O}_3$ is a straight line with increasing slope.

5. Why rate of reaction is fastest, if distilled water is not added to 50 mL (i.e D flask) 0.1 M $\text{Na}_2\text{S}_2\text{O}_3$?

Ans. Because the Concⁿ of $\text{Na}_2\text{S}_2\text{O}_3$ increases

Rate $\propto [\text{conc}]$

Remark and sign of teacher:

Observation Table:

1	Mass of 0.2 M CuSO ₄ solution	-	25 g
2	Mass of zinc metal	-	1 g
3	Total mass of solution	m	26 g
4	Initial temperature of CuSO ₄ solution	t ₁	27 °C
5	Final temperature of mixture	t ₂	37 °C
6	Rise in temperature	Δt = 37 °C (t ₂) - 27 °C (t ₁)	10 °C
7	Specific heat capacity of water	S	4.184 J/g °C

Calculation:

Heat evolved in displacement reaction is given by

$$Q = m S \Delta t$$

$$Q = 26 \times 4.184 \text{ J/g °C} \times 10 \cdot \Delta t$$

$$Q = 1087.8 \text{ J}$$

1000 mL 0.2 M CuSO₄ solution contain 0.2 mol of 'Cu'

25 mL 0.2 M CuSO₄ solution contain

$$= \frac{25 \times 0.2}{1000}$$

$$= 5 \times 10^{-3} \text{ mol of 'Cu'}$$

..... (Q) is the heat evolved, when 5×10^{-3} mol of 'Cu' is displaced it

Therefore enthalpy of displacement is

$$\Delta H = \frac{1087.8 \text{ (Q)}}{5 \times 10^{-3}}$$

$$\Delta H = 217.56 \times 10^3 \text{ J/mol}$$

$$\Delta H = 217.56 \text{ kJ/mol}$$

Space for log calculation

Number	log
--------	-----

Result:

The enthalpy of displacement of 'Cu' from copper sulphate by zinc dust $\Delta H = 217.56 \text{ kJ/mol}$
 (Negative sign is because reaction is exothermic in nature)

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. The number of water molecules present in hydrated CuSO_4 is/are

- a. 1
- b. 2
- c. 3
- d. 5

2. The metal used for connection of railway tracks after a particular distance

- a. Zn
- b. Al
- c. K
- d. Fe

3. A thermodynamic state function is a quantity

- a. used to determine heat changes
- b. whose value is independent of path
- c. used to determine pressure volume work
- d. whose value depends on temperature only

4. The enthalpies of all elements in their standard state are

- a. unity
- b. zero
- c. < 0
- d. different for each element

Short answer questions

1. Why does Zinc displace copper from CuSO_4 solution?

Ans. Zinc displaces copper from CuSO_4 solⁿ because zinc is having large value of standard oxidation potential.

[$E_{\text{oxidation}} = +0.763 \text{ V}$] than copper $E^{\circ}_{\text{oxidation}} = 0.33$

2. What is enthalpy of displacement?


Ans. The enthalpy change that occurs when one mole of substance is completely displaced from its aqueous solⁿ is called enthalpy of displacement.

3. What is specific heat of a substance?

Ans. The specific heat of substance is the amount of energy required to raise the temp of 1 gram of the substance by 1°C unit $\text{J/g}^{\circ}\text{C}$ or $\text{cal/g}^{\circ}\text{C}$.

4. Write balanced chemical equation, when Iron filings are kept in blue coloured solution of copper sulphate.

Ans. When iron filling are kept in blue coloured solⁿ of CuSO_4 , then blue colour slowly becomes colourless due to formation of FeSO_4 .



Remark and sign of teacher:

Compound No.1

Test	Observation	Inference
1. Test for Carboxylic acid group Substance + 10% NaHCO ₃ solution.	Substance soluble with brisk effervescence of CO ₂ gas.	Carboxylic group present (-C-OH)

Result: The given organic compound no.1 contains following functional group.

Name of the functional group	Structure of functional group
Carboxylic.....	-C(=O)-OH

Compound No.2

Test	Observation	Inference
1. Test for Carboxylic acid group Substance + 10% NaHCO ₃ solution.	No brisk effervescence of CO ₂ gas.	Carboxylic group absent.
2. Test for Phenolic group Substance + little water + neutral FeCl ₃ solution	green colouration.	Phenolic (Ar-OH) group present

Result: The given organic compound no.2 contains following functional group.

Name of the functional group	Structure of functional group
Phenolic.....	Ar-OH.

Compound No.3

Test	Observation	Inference
1. Test for Carboxylic acid group Substance + 10% NaHCO ₃ solution.	No brisk effervescence of CO ₂ gas.	Carboxylic group absent.
2. Test for Phenolic group Substance + little water + neutral FeCl ₃ solution	No colouration	Phenolic group absent.
3. Test for Aldehyde group (liquid) Substance + 2-3 mL Schiff's reagent/ Tollen's reagent and heat	Pink colour appears.	Aldehyde (-C=O) group present.

Result: The given organic compound no.3 contains following functional group.

Name of the functional group	Structure of functional group
Aldehyde.....	-C(=O)H

Compound No.4

Test	Observation	Inference
1. Test for Carboxylic acid group Substance + 10% NaHCO ₃ solution.	No brisk effervescence of CO ₂ gas.	Carboxylic group absent.

2. Test for Phenolic group Substance + little water + neutral FeCl_3 solution	NO colouration	phenolic group absent.
3. Test for Aldehyde group (liquid) Substance + 2-3 mL Schiff's reagent/ Tollen's reagent and heat	NO pink colouration	Aldehyde group absent
4. Test for Ketonic group (liquid) Substance + few drops of NaOH + few drops of sodium nitroprusside soln	Red colour.	Ketonic group present ($>\text{C}=\text{O}$)

Result: The given organic compound no.4 contains following functional group.

Name of the functional group	Structure of functional group
Ketonic	$>\text{C}=\text{O}$.

Compound No.5

Test	Observation	Inference
1. Test for Carboxylic acid group Substance + 10% NaHCO_3 solution.		
2. Test for Phenolic group Substance + little water + neutral FeCl_3 soln		
3. Test for Aldehyde group (liquid) Substance + 2-3 mL Schiff's reagent/ Tollen's reagent and heat		
4. Test for Ketonic group (liquid) Substance + few drops of NaOH + few drops of sodium nitroprusside soln		
5. Test for Amino group Substance + conc. HCl and shake to dissolve. Cool under tap water, add excess of NaNO_2 + β -naphthol in NaOH		

Result: The given organic compound no.5 contains following functional group.

Name of the functional group	Structure of functional group
.....

Remark and sign of teacher:

MCO

Select [✓] the most appropriate answer from given alternatives of each sub question.

- Appropriate answer from given alternatives of each one question.

 - The functional group present in Methanal is
 - a. $\text{-C}=\text{O}$
 - b. $\text{C}=\text{O}$
 - c. $\text{-C}(\text{OH})=\text{O}$
 - d. $\text{-N}(\text{H})\text{H}$
 - Identify the functional group of primary amine
 - a. -NH_2
 - b. $\text{-N}(\text{H})\text{H}$
 - c. $\text{N}(\text{H})_2$
 - d. NH_3
 - The distinguishing test between alcohol and phenol is -
 - a. bicarbonate test
 - b. neutral FeCl_3 test
 - c. sodium nitroprusside test
 - d. tollen's reagent test
 - Primary amines on treatment with excess of NaNO_2 , concentrated HCl and β -naphthol in alkaline medium gives.....
 - a. green colour
 - b. orange colour
 - c. violet colour
 - d. red colour
 - Acetaldehyde when heated with ammonical silver nitrate solution gives
 - a. red ppt
 - b. silver mirror
 - c. red colour
 - d. grey mirror

Short answer questions

1. Define: Functional group.
Ans..... Functional group is defined as an atom or group of atoms present in a molecule which represents its properties & responsible for characteristic reaction.


2. Draw structure of an organic compound having two functional group in molecule.
Ans.....



3. Give distinguishing test for carboxylic acid and phenol.
Ans..... Sodium bicarbonate NaHCO_3 readily reacts with carboxylic acid but does not react with phenol.

4. If alkaline KMnO_4 solution is added to an organic compound, it decolorises. Identify the functional group present in the compound?
Ans..... If alkaline KMnO_4 soln is added to an organic compound, it decolorises. It indicates that such compound contain unsaturation $\text{>C=C<} \text{ or } \text{-C=C-}$

5. Write a chemical reaction for Tollen's reagent test with ethanal.
Ans.....



Remark and sign of teacher:

Aim: To detect the presence of carbohydrates or proteins or fats and oils in the given food stuff.

Food stuff for carbohydrates (Perform any two tests)

Sr. No.	Test	Observation	Inference
1	Solubility: Substance + Water	Soluble in water	Glucose present
2	Fehling Test: Little substance + 1 mL Fehling solution 'A' and 'B'. Heat the test tube on a water bath	Red ppt. is formed	Glucose Present
3	Benedict's Test: Little substance + add 2 mL Benedict's reagent. Heat gently.	Red ppt is obtained	Carbohydrate is Present

Result : The given Sample food stuff contains carbohydrates

Food stuff for protein (Perform any two tests)

Sr. No.	Test	Observation	Inference
1	Biuret Test: Substance + 10% NaOH + few drops of 1% CuSO ₄ solution, shake well	Bluish violet colour	Protein is present
2	Xanthoprotic Test: Substance + few drops of conc. HNO ₃ . Shake a little and leave undisturbed for some time	Yellow colouration	Protein is present
3	Ninhydrin Test: Substance + NaOH solution + few drop of Ninhydrin reagent. Heat.	Intense blue colour	Protein is present

Result: The given Sample food stuff contains protein.

Food stuff for fats and oils (Perform any two tests)

Sr. No.	Test	Observation	Inference
1	Solubility: a. Substance + Water	immiscible	oil is present
	b. Substance + chloroform	miscible	oil is present
2	Spot Test: Put a drop of substance on filter paper	translucent spot increases on drying filter paper	oil is present
3	Acrolein Test: Substance + few crystals of potassium bisulphite. Heat	pungent irritating smell of formation of acrolein	Oil Present

Result: The given Sample food stuff contains fats or oils.

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. Fats and oils are -----
 ✓ a. triesters of glycerol and fatty acids b. diesters of glycerol and fatty acids
 c. glycerol and fatty acids d. glycol and fatty acids
2. ----- are insoluble in water but soluble in acidic or alkaline solutions
 ✓ a. fats b. oils c. proteins d. carbohydrates
3. Proteins contain ----- linkage.
 ✓ a. polypeptide b. polyester c. polypropylene d. polystyrene
4. Glucose powder when heated with Fehling solution gives-----
 a. blue ppt ✓ b. red ppt c. green ppt. d. black ppt

Short answer questions

1. What are Carbohydrates?


Ans.... Optically active polyhydroxy aldehydes or polyhydroxy ketones are called Carbohydrates.

2. Write the name of chemical test, when protein in food stuff is treated with concentrated HNO_3 .

Ans.... When protein is food stuff is treated with conc. HNO_3 it gives yellow colour or ppt. This called Xanthoprotein test

3. Write structure of peptide linkage present in proteins.

Ans.....



4. Mention Biuret test for proteins in short.

Ans.... Substance containing protein is treated with 1% NaOH & few drops of 1% CuSO_4 soln on shaking it gives

5. Write an important difference between oils and fats. Violet colour.

Ans.... At ordinary temp. oils are liquids while fats are solid. Oils are triesters of unsaturated fatty acids while fats are triesters of saturated fatty acids.

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. The anion does NOT give positive test with dil. H_2SO_4 and conc. H_2SO_4
 - Cl^-
 - CO_3^{2-}
 - CH_3COO^-
 - SO_4^{2-} ✓
2. When moisten oxalic acid is taken between the fingers, it smells like vinegar, indicates the presence of following anion
 - sulphate
 - nitrate
 - acetate ✓
 - nitrite
3. Canary yellow ppt. of ammonium phosphomolybdate indicates the presence of radical
 - Cl^-
 - SO_4^{2-}
 - PO_4^{3-} ✓
 - NO_3^-
4. Smell of rotten eggs of a substance with dil. H_2SO_4 indicates the presence of radical
 - SO_3^{2-}
 - SO_4^{2-}
 - S^{2-} ✓
 - $C_2O_4^{2-}$

Short answer questions

1. Write the formula of the compound formed in brown ring test.

Ans..... The compound formed in brown ring test is



2. Write the name of the complex formed, when sodium nitroprusside is added to O.S. of mixture in the confirmatory test of sulphide ion.

Ans..... When Sodium Nitroprusside is added to O.S. of mixture in the confirmatory test of S^{2-} ion, then Sodium pentacyanoferrate (II) Sulphide $Na_4[Fe(CN)_5S] \cdot S^{2-}$ is formed.

3. Write the names of anion detected with the help of dil. H_2SO_4 .

Ans..... CO_3^{2-} (Carbonate ion), NO_3^- Nitrite, S^{2-} (Sulphide ion)
 SO_3^{2-} (Sulphite ion)

4. Write the names of anion detected with the help of conc. H_2SO_4 .

Ans..... Cl^- (Chloride), Br^- (Bromide), I^- (Iodide), NO_3^- (Nitrate)
 Oxalate & CH_3COO^- (Acetate ion)

5. What is lime water and what happens on passing CO_2 gas through it?

Ans..... Lime water is solution of Calcium hydroxide in H_2O which is alkaline & turns milky as passing CO_2 gas through it.

6. Describe the layer test for iodide ions?

Ans..... Take soln containing iodine ions & Cl_2 water & chloroform
 Shake well. Chloroform layer acquires violet layer.

7. How do you test the presence of sulphide ion in dry test?

Ans..... Mixture of BaCl₂ is treated with dil. H_2SO_4 which evolves colourless gas with smell of rotten eggs (CH_3S) which is turns lead acetate paper black indicates presence of

Remark and sign of teacher:

Experiment No. 20 (Mixture No. 1)

Date: / /

Aim : Analyse two acidic (anion) radicals qualitatively from given inorganic mixture.

Apparatus : Test tubes, test tube holder, test tube stand, filter paper etc.

A. Preliminary tests:

Test	Observation	Inference
Color	White	Chloride, Nitrates, Sulphates, Carbonates of Ca^{2+} , Ba^{2+} , Sr^{2+} or Mg^{2+} may be present
Nature	Crystalline	Chloride, Br^- , I^- , NO_2^- , SO_4^{2-} , CO_3^{2-} , CH_3COO^- or compounds of

B. Dry Tests for Acidic radicals (Test tube must be dry for following tests) NH_4^+ or K^+ may be present

Test	Observations	Inference
1. Heating in a dry test tube:- Take a small quantity of the mixture in a clean and dry test tube and heat it strongly	Colourless, odourless gas, which turns lime water milky.	CO_3^{2-} , $\text{C}_2\text{O}_4^{2-}$ may be present
2. Action of dil. H_2SO_4 : Take a small quantity of mixture + dilute H_2SO_4	Evolution of colourless CO_2^- may be present Odourless gas with brisk effervescence of CO_2 gas	CO_3^{2-} may be present
3. Action of conc. H_2SO_4 : Take a small quantity of mixture + conc. H_2SO_4	colourless gas having Pungent smell	Cl^- may be present
4. Action of Cu foils and conc. H_2SO_4 : Mixture + Cu filings and conc. H_2SO_4 , heat it strongly	No brown fumes	NO_3^- is absent
5. Action of MnO_2 and conc. H_2SO_4 : Mixture + MnO_2 powder + conc. H_2SO_4	colourless gas	Cl^- may be present.

Individual dry tests for SO_4^{2-} & PO_4^{3-} (if required)

1. Test for sulphate Mixture + dil. HCl, boil the solution + BaCl_2 solution	No white ppt	SO_4^{2-} is absent
2. Test for phosphate Mixture + conc. HNO_3 , boil + excess of ammonium molybdate solution	No Canary Yellow ppt	PO_4^{3-} is absent.

C. Preparation of original solution (O.S.) : Mixture is dissolved in 20 mL quantity of distilled water in a beaker, stir with glass rod to dissolve the mixture. Clear solution is obtained. Use this O.S. for further tests of acidic radicals.

D. Wet Tests for Anion (Acidic Radicals)

1. O.S. + AgNO_3	White ppt insoluble in dil. HNO_3	Halides may be Cl^- , Br^- , I^- present.
If white ppt obtained in above test, perform following distinction test		
Distinction between Cl^- , Br^- and I^- : O.S. + dil. H_2SO_4 (till acidic) + chloroform + Cl_2 water (fresh) in excess, shake vigorously and observe the colour of chloroform layer carefully.	chloroform layer COLOURLESS	Cl^- Present
2. O.S. + $\text{Ba}(\text{NO}_3)_2$ solution	white ppt. Soluble in dil. HCl or H_2SO_4 with effervescence	CO_3^{2-} Present
3. O.S. + dil. acetic acid + (Freshly prepared) FeSO_4 solution	NO brown colouration	NO_3^- absent
4. O.S. + diphenyl amine + conc. H_2SO_4	NO evolution of NO_2 gas	NO_3^- absent
5. O.S. + FeCl_3 solution	No. PPT	PO_4^{3-} & CH_3COO^- absent
6. O.S. + dil. acetic acid + CaCl_2	NO. White PPT.	$\text{C}_2\text{O}_4^{2-}$ absent
7. O.S. + dil. H_2SO_4 + 2-3 drops of KMnO_4	Pink colour of KMnO_4 remains as it is.	SC_3^{2-} absent.

E. Confirmatory tests for Acidic Radicals

1. C.T. for first detected acidic radical ... CO_3^{2-}

Test	Observation	Inference
1. O.S + 2-3 drops Phenolphthalein indicator	Pink colour	CO_3^{2-} confirmed
2. O.S + dil HCl	Effervescence of CO_2 which turns lime water milky.	CO_3^{2-} confirmed

2. C.T. for second detected acidic radical ... Cl^-

Test	Observation	Inference
1. O.S + Lead acetate soln	White PPT.	Cl^- confirmed.
2. O.S + AgNO_3	White PPT Soluble in NH_4OH	Cl^- confirmed

Result :

The given inorganic mixture no. 1 contains following two anions (Acidic radicals)

- i) CO_3^{2-} (Carbonate)
- ii) Cl^- (Chloride)

Remark and sign of teacher:

Experiment No. 21 (Mixture No. 2)

Date: / /

Aim : Analyse two acidic (anion) radicals qualitatively from given inorganic mixture.

Apparatus : Test tubes, test tube holder, test tube stand, filter paper etc.

A. Preliminary tests:

Test	Observation	Inference
Color	Blue	chlorides, Nitrates, Sulphates, carbonates of Cu^{2+} , Cr^{2+} , Ni^{2+} may be present
Nature	crystalline	Cl^- , Br^- , I^- , NO_2^- , SO_4^{2-} , CO_3^{2-} CH_3COO^- or compounds of NH_4^+ or K^+ may be present.

B. Dry Tests for Acidic radicals (Test tube must be dry for following tests)

Test	Observations	Inference
1. Heating in a dry test tube:- Take a small quantity of the mixture in a clean and dry test tube and heat it strongly	Greenish Yellow gas which turns Starch Iodide Paper blue.	Cl^- may be Present
2. Action of dil. H_2SO_4 : Take a small quantity of mixture + dilute H_2SO_4	NO colour gas.	CO_3^{2-} , S^{2-} , NO_2^- SO_3^{2-} , CH_3COO^- absent
3. Action of conc. H_2SO_4 : Take a small quantity of mixture + conc. H_2SO_4	Colourless gas having Pungent smell	Cl^- may be Present
4. Action of Cu foils and conc. H_2SO_4 Mixture + Cu filings and conc. H_2SO_4 , heat it strongly	NO brown fumes	NO_3^- is absent
5. Action of MnO_2 and conc. H_2SO_4 Mixture + MnO_2 powder + conc. H_2SO_4	greenish yellow gas turning moist blue litmus paper red & then bleaches	Cl^- may be Present

Individual dry tests for SO_4^{2-} & PO_4^{3-} (if required)

1. Test for sulphate Mixture + dil. HCl, boil the solution + BaCl_2 solution	white ppt. insoluble in conc. HCl	SO_4^{2-} maybe Present
2. Test for phosphate Mixture + conc. HNO_3 , boil + excess of ammonium molybdate solution	No canary yellow ppt.	PO_4^{3-} is Absent.

C. Preparation of original solution (O.S) : Mixture is dissolved in 20 mL quantity of distilled water in a beaker, stir with glass rod to dissolve the mixture. Clear solution is obtained. Use this O.S. for further tests of acidic radicals.

D. Wet Tests for Anion (Acidic Radicals)

1. O.S. + AgNO_3	White ppt insoluble Halides may be present (Cl, Br, I) present in dil. HNO_3
If white ppt obtained in above test, perform following distinction test	
Distinction between Cl ⁻ , Br ⁻ and I ⁻ :	
O.S. + dil. H_2SO_4 (till acidic) + chloroform + Cl ₂ water (fresh) in excess, shake vigorously and observe the colour of chloroform layer carefully.	Chloroform layer colourless
2. O.S. + $\text{Ba}(\text{NO}_3)_2$ solution	White ppt insoluble SO_4^{2-} present in dil. HNO_3
3. O.S. + dil. acetic acid + (Freshly prepared) FeSO_4 solution	No brown colouration
4. O.S. + diphenyl amine + conc. H_2SO_4	No evolution of NO_2 gas
5. O.S. + FeCl_3 solution	No ppt
6. O.S. + dil. acetic acid + CaCl_2	No white ppt
7. O.S. + dil. H_2SO_4 + 2-3 drops of KMnO_4	Pink colour of KMnO_4 remains as it is.

E. Confirmatory tests for Acidic Radicals

1. C.T. for first detected acidic radical ... Cl⁻ ..

Test	Observation	Inference
1. O.S + Lead acetate soln.	White ppt.	Cl ⁻ confirmed
2. O.S + AgNO_3	White ppt Soluble in NH_4OH	Cl ⁻ confirmed

2. C.T. for second detected acidic radical ... SO_4^{2-} (Sulphate)

Test	Observation	Inference
1. O.S + Lead acetate soln.	White ppt. Insoluble in dil HNO_3 but soluble in hot ammonium acetate soln.	SO_4^{2-} confirmed.
2. O.S + BaCl_2	White Ppt of BaSO_4	SO_4^{2-} confirmed

Result :

The given inorganic mixture no.2 contains following two anions (Acidic radicals)
 i) Cl⁻ (Chloride)
 ii) SO_4^{2-} (Sulphate)

Remark and sign of teacher:

Experiment No. 22 (Mixture No. 3)

Date: / /

Aim : Analyse two acidic (anion) radicals qualitatively from given inorganic mixture.

Apparatus : Test tubes, test tube holder, test tube stand, filter paper etc.

A. Preliminary tests:

Test	Observation	Inference
Color	White	Chlorides, Nitrates, Sulphates, Carbonates of Ca^{2+} , Ba^{2+} , Sr^{2+} , or Mg^{2+} may be present
Nature	crystalline	Cl^- , Br^- , I^- , NO_2^- , SO_4^{2-} , CO_3^{2-} , CH_3COO^- or compounds of NH_4^+ or K^+ may be present.

B. Dry Tests for Acidic radicals (Test tube must be dry for following tests)

Test	Observations	Inference
1. Heating in a dry test tube:- Take a small quantity of the mixture in a clean and dry test tube and heat it strongly	Colourless, odourless CO_3^{2-} , $\text{C}_2\text{O}_4^{2-}$ gas which turns lime water milky.	CO_3^{2-} may be present.
2. Action of dil. H_2SO_4 : Take a small quantity of mixture + dilute H_2SO_4	Evolution of colourless, odourless gas with brisk effervescence.	CO_3^{2-} may be present.
3. Action of conc. H_2SO_4 : Take a small quantity of mixture + conc. H_2SO_4	Yellowish brown gas having pungent smell turns starch paper yellow.	Br^- may be present.
4. Action of Cu foils and conc. H_2SO_4 : Mixture + Cu filings and conc. H_2SO_4 , heat it strongly	NO Brown fumes	NO_3^- is absent
5. Action of MnO_2 and conc. H_2SO_4 : Mixture + MnO_2 powder + conc. H_2SO_4	Reddish brown gas turning starch paper yellowish brown	Br^- may be present.

Individual dry tests for SO_4^{2-} & PO_4^{3-} (if required)

1. Test for sulphate Mixture + dil. HCl, boil the solution + BaCl_2 solution	NO white PPT	SO_4^{2-} is absent.
2. Test for phosphate Mixture + conc. HNO_3 , boil + excess of ammonium molybdate solution	No Canary Yellow PPT.	PO_4^{3-} is absent.

C. Preparation of original solution (O.S) : Mixture is dissolved in 20 mL quantity of distilled water in a beaker, stir with glass rod to dissolve the mixture. Clear solution is obtained. Use this O.S. for further tests of acidic radicals.

D. Wet Tests for Anion (Acidic Radicals)

1. O.S. + AgNO_3 ,	Pale yellow ppt Partly soluble in NH_4OH	Br^- may be present
-----------------------------	---	------------------------------

If white ppt obtained in above test, perform following distinction test

Distinction between Cl^- , Br^- and F^- :	Chloroform layer Yellowish brown.	Br^- present
O.S. + dil. H_2SO_4 (till acidic) + chloroform + Cl_2 water (fresh) in excess, shake vigorously and observe the colour of chloroform layer carefully.	White ppt Soluble in dil HCl	CO_3^{2-} present.
3.O.S. + dil. acetic acid + (Freshly prepared) FeSO_4 solution	NO brown colouration	NO_2^- absent.
4. O.S. + diphenyl amine + conc. H_2SO_4	NO evolution of NO_2 gas.	NO_3^- absent
5. O.S. + FeCl_3 solution	NO ppt	PO_4^{3-} & CH_3COO^- absent
6. O.S. + dil. acetic acid + CaCl_2	NO white ppt.	$\text{C}_2\text{O}_4^{2-}$ absent
7. O.S. + dil. H_2SO_4 + 2-3 drops of KMnO_4	Pink colour of	SO_4^{2-} absent.

E. Confirmatory tests for Acidic Radicals

1. C.T. For first detected acidic radical ... CO_3^{2-} same as it

Test	Observation	Inference
1. O.S + 2-3 drops of Phenolphthalein indicator.	Pink colour	CO_3^{2-} confirmed
2. O.S + dil. HCl.	Effervescence of CO_2 which turns lime water milky	CO_3^{2-} confirmed.

2. C.T. For second detected acidic radical ... Br^- (Bromide)

Test	Observation	Inference
1. O.S + MnO_2 + conc. H_2SO_4 Heat.	Raddish brown gas turning Starch paper Yellowish brown.	Br^- confirmed.
2. O.S + Cl_2 water + CCl_4	Chloroform layer aquires Yellow colour.	Br^- confirmed.

Result :

The given inorganic mixture no. 3 contains following two anions (Acidic radicals)
 i) CO_3^{2-} (Carbonate)
 ii) Br^- (Bromide)

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. A colourless solid 'A' produces black spots on the skin. It's aqueous solution gives brown ring test and also gives yellow ppt. with KI solution 'A' could be
 - a. copper nitrate
 - b. zinc nitrate
 - c. silver nitrate
 - d. lead nitrate
2. Which of the following sulphide is completely precipitated in acidic medium
 - a. HgS
 - b. PbS
 - c. CdS
 - d. CuS
3. The colour of precipitate observed, when potassium iodide is added to water soluble lead salts is
 - a. black
 - b. white
 - c. yellow
 - d. red
4. Which of the following is NOT a preliminary test
 - a. charcoal cavity test
 - b. flame test
 - c. NaOH test
 - d. brown ring test
5. The colour of carbonate cation precipitate of group V is.....
 - a. black
 - b. green
 - c. white
 - d. yellow

Short answer questions

1. The removal of H_2S gas from filtrate is must before doing analysis of group III and V.

Why?

Ans. If H_2S gas isn't boiled off then group IV cations like Cu^{2+} get precipitated. It would react with HNO_3 & would be oxidised to colloidal Sulphur which would interfere with further analysis.

2. Write the name of brown ppt. obtained, when Nessler's reagent is added in O.S. for confirmatory test of ammonium ion.

Ans. Brown Ppt. = Basic mercury (II) amido-iodine

3. During the detection of groups, if Group-I is detected then for next Group test what care should be taken?

Ans. If group-I is detected then reagent must be added in bulk O.S. for complete ppt formation of group I cation then filter it & filtrate is taken for further tests.

4. Why cracking noise/decrepitation is observed during dry test in qualitative analysis of cations?

Ans. Lead Nitrate, Barium Nitrate, Potassium bromide, Sodium chloride etc. make cracking noise due to decomposition on heating in a dry test tube. Also due to heating bigger crystal breaks up into smaller ones.

5. Name a cation, which is not obtained from metal.

Ans. cation Ammonium ion NH_4^+ is not obtained from metal

Remark and sign of teacher:

Experiment No. 23 (Mixture No. 1)

Date: / /

Aim : Analyse two basic (cation) radicals qualitatively from given inorganic mixture.

Apparatus : Test tubes, test tube holder, test tube stand, filter paper etc.

A. Preliminary tests:

Test	Observation	Inference
Color	Blue	Cu^{2+} may be Present
Nature	Crystalline	Water soluble salts like Nitrates, Sulphates & halides

B. Dry Tests for Basic radicals (Test tube must be dry for this test)

Test	Observation	Inference
1. Heating in a dry test tube:- Take a small quantity of the mixture a clean and in dry test tube and heat it strongly.	Formation of White Sublimate	NH_4^+ may be Present
2. Charcoal Cavity Test: Mixture + Na_2CO_3 solid in 1:2 proportion placed in fresh charcoal cavity, moisten with a drop of water. Heat it with blow pipe in a reducing (yellow) flame	Substance fuses & sinks in the cavity	NH_4^+ or Ca^{2+} salts may be present
3. NaOH Test: Mixture + NaOH solution & heat. Hold moist turmeric paper near the mouth of the test tube	Moist turmeric Paper turns brown.	NH_4^+ may be Present.
4. Flame Test: Prepare a paste of the given mixture with conc. HCl on a watch glass. Make a small loop at the end of the platinum wire & dip it in the mixture or use glass rod. Heat it on oxidising flame (Blue) observe the colour change of the flame.	Bluish white	Cu^{2+} may be Present

C. Preparation of original solution (O.S.)

Take a small quantity of mixture in a beaker add 20 mL of distilled water, stir with glass rod to dissolve the mixture. If mixture does not dissolve completely then warm it to dissolve. Clear solution is obtained, which is used as a O.S for further tests.

1. Analysis of Group zero (NH_4^+)

Test	Observation	Inference
1. O.S. + NaOH solution + Heat, test with moist turmeric paper.	Evolution of NH_3 gas which turns moist turmeric paper brown.	Group zero Present (NH_4^+)

2. O.S. + NaOH solution + Heat, Bring a glass rod dipped in conc. HCl. near the mouth of the test tube.	Dense white fumes of NH_4Cl .	NH_4^+ Present.
C.T. for NH_4^+ : i. O.S. + Nessler's reagent in excess	Brown ppt	NH_4^+ confirmed.
ii. O.S + Picric Acid (2, 4, 6 trinitro phenol)	Yellow ppt of ammonium picrate	NH_4^+ confirmed.

Detection of group (if group zero is absent, two groups must be detected following test)

Test	Observation	Inference
1. O.S. + dil. HCl	NO ppt.	Group I Absent
2. O.S./Filtrate + dil HCl (heat) + H_2S gas or water	Black ppt.	Group II ($\text{Cu}^{2+}, \text{As}^{3+}$) present
3. O.S./Filtrate (Remove H_2S) + NH_4Cl (equal) + NH_4OH (till alkaline to litmus)	NO ppt	Group III Absent
4. O.S./Filtrate + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + H_2S gas or water	NO ppt	Group IV Absent
5. O.S./Filtrate (Remove H_2S) + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + $(\text{NH}_4)_2\text{CO}_3$	NO ppt	Group V absent
6. O.S./Filtrate + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + Na_2HPO_4	NO ppt	Group VI absent

2. Analysis of first detected group

Analysis of II group.

Test	Observations	Inference
O.S + dil. HCl + H_2S gas or water.	Black ppt	Pb^{2+} or Cu^{2+} may be present.
Group II Black ppt H:1 HNO_3 Ppt dissolves + dil. H_2SO_4	NO ppt.	Cu^{2+} present.
Above Soln + NH_4OH in excess,	Blue colour.	Cu^{2+} present.

3. C.T. for first detected radical C.T. For Cu^{2+} .

Test	Observations	Inference
1. Above soln + acetic acid + $\text{K}_4\text{Fe}(\text{CN})_6$	chocolate brown ppt	Cu^{2+} confirmed.
2. Above soln + KI.	White ppt in brown coloured solution.	Cu^{2+} confirmed.

2. Analysis of second detected group

Test	Observations	Inference
1.		
2.		
3.		

3. C.T. for second detected radical

Test	Observations	Inference
1.		
2.		

Result:-

The given inorganic mixture no. 1 contains following two cations (Basic Radicals)

- i) Cu^{2+} (Copper)
- ii) NH_4^+ (Ammonium)

Remark and sign of teacher:

Experiment No. 24 (Mixture No. 2)

Date: / /

Aim : Analyse two basic (cation) radicals qualitatively from given inorganic mixture.

Apparatus : Test tubes, test tube holder, test tube stand, filter paper etc.

A. Preliminary tests:

Test	Observation	Inference
Color	Green.	Ni^{2+} may be Present
Nature	crystalline	Water soluble Salts like Nitrates

B. Dry Tests for Basic radicals (Test tube must be dry for this test)

Test	Observation	Inference
1. Heating in a dry test tube:- Take a small quantity of the mixture a clean and in dry test tube and heat it strongly.	Formation of white sublimate	NH_4^+ may be present
2. Charcoal Cavity Test: Mixture + Na_2CO_3 solid in 1:2 proportion placed in fresh charcoal cavity, moisten with a drop of water. Heat it with blow pipe in a reducing (yellow) flame	Substance fuses & sinks in the cavity	NH_4^+ or Cat^{2+} Salts may be present.
3. NaOH Test: Mixture + NaOH solution & heat. Hold moist turmeric paper near the mouth of the test tube	Moist turmeric Paper turns brown.	NH_4^+ may be present
4. Flame Test: Prepare a paste of the given mixture with conc. HCl on a watch glass. Make a small loop at the end of the platinum wire & dip it in the mixture or use glass rod. Heat it on oxidising flame (Blue) observe the colour change of the flame.	Bluish green	Zn^{2+} & Mn^{2+} Salts may be present.

C. Preparation of original solution (O.S.)

Take a small quantity of mixture in a beaker add 20 mL of distilled water, stir with glass rod to dissolve the mixture. If mixture does not dissolve completely then warm it to dissolve. Clear solution is obtained, which is used as a O.S for further tests.

1. Analysis of Group zero (NH_4^+)

Test	Observation	Inference
1. O.S. + NaOH solution + Heat, test with moist turmeric paper.	Evolution of NH_3 gas which turns moist turmeric Paper 114 brown	Group zero present (NH_4^+)

2. O.S. + NaOH solution + Heat, Bring a glass rod dipped in conc. HCl. near the mouth of the test tube.	Dense white fumes of NH_4Cl	NH_4^+ Present
C.T. for NH_4^+ : i. O.S. + Nessler's reagent in excess	Brown ppt	NH_4^+ confirmed.
ii. O.S. + Picric Acid (2, 4, 6 trinitro phenol)	Yellow ppt of ammonium picrate	NH_4^+ confirmed.
Detection of group (if group zero is absent, two groups must be detected following test)		
Test	Observation	Inference
1. O.S. + dil. HCl	NO white ppt	Group I absent
2. O.S./Filtrate + dil HCl (heat) + H_2S gas or water	NO Black ppt	Group II absent
3. O.S./Filtrate (Remove H_2S) + NH_4Cl (equal) + NH_4OH (till alkaline to litmus)	NO ppt	Group III absent
4. O.S./Filtrate + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + H_2S gas or water	Black ppt	Group IV Present
5. O.S./Filtrate (Remove H_2S) + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + $(\text{NH}_4)_2\text{CO}_3$	NO ppt	Group V absent
6. O.S./Filtrate + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + Na_2HPO_4	NO ppt	Group VI absent

2. Analysis of first detected group Analysis of Group IV

Test	Observations	Inference
O.S. + NH_4Cl + NH_4OH + H_2S gas	Black ppt	$\text{Ni}^{2+}, \text{CO}^{2+}$ Present
Colour of original soln	Greenish	Ni^{2+} Present.
Dissolve the Black ppt of IV in aqua regia (conc. $\text{HCl} + \text{HNO}_3$) 3:1 proportion dilute with water. Use this soln for C.T. Ni^{2+} .		

3. C.T. for first detected radical C.T. for Ni^{2+} .

Test	Observations	Inference
1. Above soln + NH_4OH + dimethyl glyoxime	Scarlet Red ppt	Ni^{2+} confirmed.
2. Above soln + NaOH	Light green ppt	Ni^{2+} confirmed.

2. Analysis of second detected group

Test	Observations	Inference
1.		
2.		
3.		

3. C.T. for second detected radical

Test	Observations	Inference
1.		
2.		

Result:-

The given inorganic mixture no. 2 contains following two cations (Basic Radicals)

- i) Ni^{2+} (Nickelion)
- ii) NH_4^+ (Ammonium)

Remark and sign of teacher:

Experiment No. 25 (Mixture No. 3)

Date: / /

Aim : Analyse two basic (cation) radicals qualitatively from given inorganic mixture.

Apparatus : Test tubes, test tube holder, test tube stand, filter paper etc.

A. Preliminary tests:

Test	Observation	Inference
Color	White	$\text{Sr}^{2+}, \text{Ba}^{2+}$ may be present
Nature	Crystalline	Water soluble salts like Nitrates, Sulphates & halides of NH_4^+ , Al^{3+} etc may be present

B. Dry Tests for Basic radicals (Test tube must be dry for this test)

Test	Observation	Inference
1. Heating in a dry test tube:- Take a small quantity of the mixture a clean and in dry test tube and heat it strongly.	Formation of white sublimate	NH_4^+ may be present
2. Charcoal Cavity Test: Mixture + Na_2CO_3 solid in 1:2 proportion placed in fresh charcoal cavity, moisten with a drop of water. Heat it with blow pipe in a reducing (yellow) flame	Substance fuses & sinks in the cavity	NH_4^+ or Ca^{2+} salt may be present
3. NaOH Test: Mixture + NaOH solution & heat. Hold moist turmeric paper near the mouth of the test tube	Moist turmeric Paper turns brown	NH_4^+ may be present
4. Flame Test: Prepare a paste of the given mixture with conc. HCl on a watch glass. Make a small loop at the end of the platinum wire & dip it in the mixture or use glass rod. Heat it on oxidising flame (Blue) observe the colour change of the flame.	Bluish green	Zn^{2+} & Mn^{2+} salt may be present

C. Preparation of original solution (O.S.)

Take a small quantity of mixture in a beaker add 20 mL of distilled water, stir with glass rod to dissolve the mixture. If mixture does not dissolve completely then warm it to dissolve. Clear solution is obtained, which is used as a O.S for further tests.

1. Analysis of Group zero (NH_4^+)

Test	Observation	Inference
1. O.S. + NaOH solution + Heat, test with moist turmeric paper.	Evolution of NH_3 gas which turns moist turmeric Paper [117] brown	Group zero present (NH_4^+)

2. O.S. + NaOH solution + Heat, Bring a glass rod dipped in conc. HCl. near the mouth of the test tube.	Dense white fumes of NH_4Cl	NH_4^+ Present
C.T. for NH_4^+ : i. O.S. + Nessler's reagent in excess	Brown Ppt	NH_4^+ confirmed
ii. O.S + Picric Acid (2, 4, 6 trinitro phenol)	Yellow ppt of ammonium picrate	NH_4^+ confirmed

Detection of group (if group zero is absent, two groups must be detected following test)

Test	Observation	Inference
1. O.S. + dil. HCl	NO Ppt	Group I st absent
2. O.S./Filtrate + dil HCl (heat) + H_2S gas or water	NO Black ppt	Group II nd absent
3. O.S./Filtrate (Remove H_2S) + NH_4Cl (equal) + NH_4OH (till alkaline to litmus)	NO Ppt.	Group III absent
4. O.S./Filtrate + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + H_2S gas or water	white ppt	Group IV Present
5. O.S./Filtrate (Remove H_2S) + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + $(\text{NH}_4)_2\text{CO}_3$	NO Ppt	Group V absent
6. O.S./Filtrate + NH_4Cl (equal) + NH_4OH (till alkaline to litmus) + Na_2HPO_4	NO Ppt	Group VI absent

2. Analysis of first detected group (IV group)

Test	Observations	Inference
O.S + $\text{NH}_4\text{Cl} + \text{NH}_4\text{OH} + \text{H}_2\text{S}$ gas	White Ppt	Zn^{2+} present.
Dissolve the white PPT in dil HCl & remove H_2S gas by boiling. Use this soln for C.T of Zn^{2+} .		

3. C.T. for first detected radical C.T. for Zn^{2+}

Test	Observations	Inference
1. O.S + NaOH	white ppt soluble in excess of NaOH but reprecipitated by H_2S	Zn^{2+} confirmed.
2. O.S + $K_4[Fe(CN)_6]$	white ppt	Zn^{2+} confirmed.

2. Analysis of second detected group

Test	Observations	Inference
1.		
2.		
3.		

3. C.T. for second detected radical

Test	Observations	Inference
1.		
2.		

Result:-

The given inorganic mixture no. 3 contains following two cations (Basic Radicals)

- i) Zn^{2+} (Zinc cation)
- ii) NH_4^+ (Ammonium)

Remark and sign of teacher:

Activity No.5

Date: / /

Aim: Identify Resin Identification Code (RIC) / Recycle Symbols of plastics on different plastic materials.

Theory: The Resin Identification Code (RIC) was developed in 1988 by the plastics Industry association. It was created for workers in the plastic and recycling industry to be able to sort and recycle plastics more efficiently. Each RIC corresponds to a specific type of resin used in a plastic product.

Chronological order of when that plastic become recyclable :

- '1' signifies that the product is made out of **Polyethylene terephthalate** (PET) (beverage bottles, cups, water bottles, soft drink, etc.)
- '2' signifies **high-density polyethylene** (HDPE) For example (shampoo bottles, water bottles, milk bottles and jugs, etc.)
- '3' signifies **Polyvinyl chloride** (PVC) For example (PVC pipes, vinyl siding, flooring, window frames, etc.)
- '4' signifies **low-density polyethylene** (LDPE) For example (plastic bags, shopping bags, soft bottles, tubing, squeezable bottles, etc.)
- '5' signifies **Polypropylene** (PP) For example (auto parts, industrial fibre, food containers, microwaveable wares, etc.)
- '6' signifies **polystyrene** (PS) For example (plastic utensils, CD cases, styrofoam, toys, hard packings, cafeteria trays etc.)
- '7' signifies other plastics, such as **acrylic, nylon, polycarbonate** and **polylactic acid** (PLA) for example. Baby feeding bottle, CD production etc.

* As per **United State Army Garrison-Hawaii** (USAG-HI) and **United State Environmental Protection Agency** (USEPA) the plastic material having number 6, 7 or with no number are not recyclable.

Apparatus: Various types of plastic materials, like mineral water bottle, transparent cell phone's back cover, PVC pipe, soft drink bottles, plastic milk bags, etc.

Procedure:

1. Take one by one given plastic material.
2. Observe the bottom or lower side or all sides of plastic material and find the triangle  with number.
3. Mention your remark in observation table whether material is recyclable or not

Observation and Result Table:

Sr. No.	Name of Plastic Materials	RIC Code	Remark (Recyclable or Not)
1	Mineral Water bottle	1	Recyclable
2	Cell phone's transparent plastic back cover	6	Not Recyclable
3	PVC pipe	3	Recyclable
4	Plastic beaker	6	Not Recyclable
5	Baby's feeding bottle	7	Not Recyclable
6	Soft drink bottle	1	Recyclable
7	Plastic container (empty)	5	Recyclable

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. The following symbol  indicate

- a. universal recycling b. universal decomposition
c. universal reuse d. use and throw

2. The RIC code on cafeteria tray is

- a. 1 b. 3 c. 4 d. 6

3. If RIC code on plastic material is 7 then which type of material will be there

- a. polylactic acid b. polyethylene terephthalate
c. high density polyethylene d. low density polyethylene

4. The RIC on auto vehicle part is

- a. 2 b. 4 c. 5 d. 7

5. The code used for recycle symbols on different plastic material is

- a. RIC b. MIC c. TLC d. PPC



Short answer questions

1. What it indicates if there is no  such symbol on plastic material?

Ans.....The plastic material are not recyclable.

2. Write RIC for beverage bottles?

Ans.....The RIC for beverage bottle is 1.

3.  signifies which type of plastic material?

Ans.....Polypropylene

4. Write RIC for food containers.

Ans.....RIC for food container is 5

5. What was the purpose for development of RIC in 1988?

Ans.....It was created for workers in Plastic & recycling industry to be able for plastic more efficiency.

Remark and sign of teacher:

Activity No.6

Date: / /

Aim: To prepare 100 mL of 0.1 M standard solution of ferrous ammonium sulphate (Mohr's salt)

Theory: The molecular formula of ferrous ammonium sulphate is $\text{FeSO}_4 \cdot (\text{NH}_4)_2\text{SO}_4 \cdot 6\text{H}_2\text{O}$.

Its molar mass is 392 g/mol.

1000 mL of 1 M solution of F.A.S. solution contains 392 g of $\text{FeSO}_4 \cdot (\text{NH}_4)_2\text{SO}_4 \cdot 6\text{H}_2\text{O}$

Hence 100 mL of 0.1 M solution of F.A.S. will contain = $\frac{100 \times 0.1 \times 392}{1000 \times 1} = 3.92 \text{ g of F.A.S.}$

Apparatus: 100 mL volumetric flask, beaker, glass rod, balance watch glass, etc.

Chemicals: F.A.S. (Mohr's salt), distilled water, dilute H_2SO_4

Procedure:

1. Weigh accurately 3.92 g Ferrous ammonium sulphate on watch glass.
2. Transfer the weighed F.A.S. to a beaker and wash the watch glass with distilled water and transfer washing to the beaker.
3. Add little distilled water to dissolve it by stirring.
4. Transfer the solution of F.A.S. from beaker to 100 mL standard flask.
5. Wash the beaker twice with distilled water and transfer washing to the standard flask.
6. Dilute the solution up to the mark on standard flask to make volume 100 mL.
7. Add few drops of H_2SO_4 to get clear solution.

Result/Conclusion: 0.1 M Standard soln of Ferrous Ammonium Sulphate is prepared.

Remark and sign of teacher:

MCQ

Select [✓] the most appropriate answer from given alternatives of each sub question.

1. The molar mass of F.A.S. is-----
 a. 392 u b. 392 g/mol c. 39.2 u d. 39.2 g/mol
2. To prepare standard 1 M solution of F.A.S.....
 a. 0.392 g of F.A.S. is dissolved in distilled water and diluted to 1 L of solution
 b. 3.92 g of F.A.S. is dissolved in distilled water and diluted to 1 L of solution
 c. 39.2 g of F.A.S. is dissolved in distilled water and diluted to 1 L of solution
 d. 392 g of F.A.S. is dissolved in distilled water and diluted to 1 L of solution
3. The second name for ferrous ammonium sulphate (F.A.S.) is
 a. Rochelle salt b. Mohr's salt c. common salt d. sea salt
4. To prepare 100 mL of standard 0.1 M solution of F.A.S.
 a. 0.392 g of F.A.S. is dissolved in distilled water and diluted to 100 cm^3
 b. 3.92 g of F.A.S. is dissolved in distilled water and diluted to 100 cm^3
 c. 39.2 g of F.A.S. is dissolved in distilled water and diluted to 100 cm^3
 d. 392 g of F.A.S. is dissolved in distilled water and diluted to 100 cm^3
5. The green colour of 0.1 M F.A.S. solution turns to which colour if we keep that solution open in an air for a long time?
 a. red b. blue c. yellowish brown d. violet

Short answer questions

1. Calculate the amount of F.A.S. required to prepare 1000 mL of 0.1 M standard solution of F.A.S..

Ans. 39.2 gm of F.A.S. is required to prepare 100 mL of 1 M soln. Hence 39.2 gm of F.A.S. is required to prepare 1000 mL of 0.1 M F.A.S. soln.

2. Calculate the molar mass of F.A.S. ($\text{Fe} = 56, \text{S} = 32, \text{N} = 14, \text{O} = 16, \text{H} = 1$)

Ans. $\text{Fe}(\text{NH}_4)_2\text{SO}_4 \cdot 6\text{H}_2\text{O}$
 $(56 + 32 + 64) + (14 + 4) \times 2 + 32 + 16 \times 4 + 12 + 96$
 $= 392 \text{ g/mol}$

3. Why the 0.1 M F.A.S. solution becomes yellowish brown on keeping open for a long time?

Ans. The 0.1 M F.A.S. soln becomes yellowish brown colour on keeping open for a long time because Fe^{2+} (green) ions undergoes oxidation to Fe^{3+} (yellow brownions).

4. Why it is necessary to add few drops of conc. H_2SO_4 in the preparation of F.A.S. solution?

Ans. It is necessary to add few drops of conc. H_2SO_4 in the preparation of F.A.S. soln because of to get the clear soln & to prevent hydrolysis of Fe^{2+} ions or to prevent

5. How many water of crystallisation are present in Mohr's salt?

Ans. In Mohr's salt 6 number of H_2O molecule are present.

Remark and sign of teacher: